




TOP NEUROPOLE Needles

For Safe Thermolesion Procedures
and Regional Anesthesia



TOP NEUROPOLE Needles

TOP Has Been Delivering These High-Quality Needles to the World for More Than 30 Years.


TOP needles are the gold standard for safe and precise Radiofrequency Lesioning (Thermolesions) and Percutaneous Nerve Blocks by injection of a local anesthetic (Regional Anesthesia).

More than 30 years of worldwide experience with Teflon-insulated TOP needles, backed by many international publications, helps to guarantee safe, comfortable and reliable procedures.

TOP NEUROPOLE Needles for Regional Anesthesia



TOP NEUROPOLE Needles for Pain Management



TOP NEUROPOLE Needles for Regional Anesthesia Ordering Information

Product Name	Product number	Gauge	Needle Length (mm) *Teflon coating	Active Tip (mm)	Package
TOP NEUROPOLE Needle ST (standard)	627081	23G	30	1 (Cutting Edge)	10 pcs./box
	627082		50		
	627083		60		
	627084		100		
	627085		150		
TOP NEUROPOLE Needle X (same wire and tubing configuration as XE Needle)	627079	23G	60	1 (Cutting Edge)	600 pcs./carton
	627080		100		

TOP NEUROPOLE Needles for Pain Management Ordering Information

Product Name	Product number	Gauge	Needle Length (mm) *Teflon coating	Active Tip (mm)	Package
TOP NEUROPOLE Needle TL (temperature control)	627066	20G	145	2	480 pcs./carton
	627067		145	5	
	627068		145	10	
	627069		145	15	
	627070	21G	100	5	10 pcs./box
	627071	22G	54	4	
	627072		100.5	2	
	627073		100.5	5	
627074	100.5		10		
TOP NEUROPOLE Needle XE (no exact temperature control)	627075	23G	60	5	600 pcs./carton
	627076		100	5	480 pcs./carton
	627077		150	5	400 pcs./carton
TOP NEUROPOLE Needle RC (radio contrast)	627078	24G	200	5	600 pcs./carton
	627086		60 (without coating)	-	
	627087	22G	100 (without coating)	-	

Note: TOP, NEUROPOLE and POLE are Trademarks by TOP Corporation.

Distributor:



TOP Corporation

19-10, Senjunakai-cho, Adachi-ku, Tokyo 120-0035, Japan
Tel: 81(Japan)-3-3882-3101 Fax: 81(Japan)-3-3881-8163

EQUIP MEDIKEY B.V.

Edisonstraat 16j - 2809 PB Gouda, The Netherlands
TEL: (0182) 573 293 Int.+31 182 573 293
FAX: (0182) 531 364 Int.+31 182 531 346
E-mail: info@equip.nl