



the key to medical innovation

Redeveloped for Pulsed Radio Frequency, with  
5 mm, 10 mm and 15 mm tip lengths

Comfortable handling and easy to maneuver

Long flexible fluid line with Luer Lock

Solid and safe electrical connection

Prepared for a new generation of PRF generators



## Neuropole XE-PRF

Needles for pulsed radio frequency  
procedures

## Top Neuropole XE-PRF Needle

The classic, all familiar XE-needles, mainly used for non-temperature-controlled lesions on the medial branch, redeveloped! It resulted in a very sophisticated new needle, specifically intended for Pulsed Radio Frequency procedures: the Neuropole XE-PRF.

With a new and improved ergonomic hub, the stiff shaft, coated in our high-grade medical Teflon and a sharper bevel, handling and maneuvering are very comfortable. The long and flexible fluid line is equipped with a Luer Lock fitting and secured cap. The electrical wire with female connector grants a safe connection to the (P)RF generator.

Like XE, the Neuropole XE-PRF needles are available in a variety of lengths, depending on the targeted nerves. New tip lengths of 10 mm and 15 mm extend that range, giving access to an entirely new set of treatment options.



Top Neuropole XE-PRF Needles - Ordering information				
Order number	Needle length	Gauge	Active tip	Package
654101	60 mm	23 G	5 mm	10 pcs/box
654104	60 mm		10 mm	
654102	100 mm		5 mm	
654105	100 mm		10 mm	
654103	150 mm		5 mm	
654106	150 mm		10 mm	
E63283	150 mm	20 G	15 mm	

