



- ▶ 150 cm unobstructed artefact free imaging area
- ▶ Smart & simple design
- ▶ Perfectly manoeuvrable
- ▶ Easy to use

## C-Table

The new C-Table has been designed for a wide range of pain management procedures. Surprisingly smart & simple.



the key to medical innovation


The new C-Table, which combines 150cm of unobstructed artefact free imaging platform, provides easy access for the most challenged patients.

Highest height: 98 cm



Lowered height: 58 cm



Carbon fibre top:  
208 cm x 61 cm



Fully locking  
medical grade  
anti-static castors

Narrow base frame for  
ease of patient transfer

### Measurements



# C-Table

The new C-Table, which combines 150cm of unobstructed artefact free imaging platform, provides easy access for the most challenged patients.

## Standard features

- ▶ 150 cm unobstructed artefact free imaging area
- ▶ Large carbon top fibre: 208 cm x 61 cm
- ▶ Lowered height: 58 cm providing easy access
- ▶ Highest height: 98 cm
- ▶ 200 kg safe working load
- ▶ Narrow base frame for ease of patient transfer
- ▶ Fully locking medical grade anti-static castors
- ▶ Hand switch
- ▶ Top frame potential grounding point
- ▶ Easy access emergency stop & isolation button for additional patient safety
- ▶ Manufactured in the UK and fully CE marked

## Perfectly manoeuvrable

- ▶ one step locking
- ▶ 360 degrees of movement
- ▶ in-line locking

## 200 kg

Safe working load of 200 kg

## STOP

Easy access emergency stop & isolation button for additional patient safety

## Available options

- ▶ Foot switch
- ▶ Pushing & positioning handrail
- ▶ Battery back up (Lithium)

# The new C-Table, surprisingly smart & simple!

<b>Codes</b>	<b>Description</b>
CTABLE001	C-Table Imaging Electric table for pain management, carbon fibre top, adjustable from 58 cm up to 98 cm

## **Optional accessories**

CTABLEAC01	Stainless steel Euro accessory rails (pair of)
CTABLEAC02	Lithium battery operating system
CTABLEAC03	Handle bar (stainless steel)
CTABLEAC05	Foot switch controls
CTABLESP08	Mattress with face hole (visco foam/anti-static)

## **Spare parts**

CTABLESP01	Carbon fibre top
CTABLESP02	LC3 400 mm stroke column
CTABLESP03	CO71 control box
CTABLESP04	Euro powerlead
CTABLESP05	HB30 handset 2 button
CTABLESP06	Foot switch controls
CTABLESP07	Lithium battery operating system
CTABLESP09	Emergency stop & isolation button
CTABLESP10	Full-length accessory rails (pair of)
CTABLESP11	Handle bar (stainless steel)
CTABLESP12	100 mm front kick-lock castors
CTABLESP13	125 mm rear linked castors
CTABLESP14	Linked castor double foot pedals LH/RH

## **More information?**

Please send us an email at [info@equip.nl](mailto:info@equip.nl)  
and we will get back to you!