

REGO'S LOWER BODY SHIELDS

X-ray protection for C-Table



High-level protection for physicians and medical personnel

- Custom-made sets of single lead strips for C-Arm imaging and treatment tables
- Easy to apply and can be mounted on both sides of the table
- Flexible, light and contain no ferromagnetic materials
- Protection level of Pb 0,50 mm
- Rich, dark blue upholstery

Safe Fluoroscopy Procedures with C-Table, *Smart & Simple*

C-Table is a *Smart & Simple* C-Arm imaging and treatment table, developed for a wide range of pain management procedures. For safe fluoroscopy procedures, we have a custom-made set of Lower Body Shields available, produced by Rego X-Ray GmbH. It offers high-level X-Ray protection with a lead equivalent value of Pb 0,50 mm.

The set is comprised of 14 lead strips, to apply on both sides of the table for maximum protection. Mounting — and rearranging, if needed — are very easy with the durable loop fastener.

On long days of diagnostics and guided needle approaches, you can rely on the shelter these well-crafted shields provide. Comforting knowledge for the full team of pain physicians and medical personnel on the floor.

Our Lower Body Shields for C-Arm Imaging Tables, standard set specifications

Included are:

Reference	Description	Qty
EQ001-75	Custom-made single lead strip Pb 0,5 mm, 75 x 19 cm, w/ loop fastener 38 x 13 cm	12 pcs
EQ001-40	Custom-made single lead strip Pb 0,5 mm, 40 x 19 cm, w/ loop fastener 38 x 13 cm	2 pcs

